

SUPPLEMENTARY SHEET - BULLS BY LOT NUMBER



Thursday 7th September 2023

LOT	TAG	SEMEN TEST	WEIGHT (KG)	BULL PEN
1	T198	PASS	452	14
2	T204	PASS	468	14
3	T182	PASS	494	14
4	T158	PASS	436	14
5	T285	PASS	495	14
6	T524	PASS	408	14
7	T162	PASS	551	14
8	T284	PASS	458	14
9	T518	PASS	459	19
10	T157	PASS	448	19
11	T236	PASS	483	19
12	T211	PASS	459	19
13	T229	WITHDRAWN	442	WITHDRAWN
14	T203	PASS	475	19
15	T163	PASS	467	19
16	S326	PASS	607	9
17	S227	PASS	521	9
18	S221	PASS	525	9
19	S420	PASS	555	9
20	S421	PASS	537	9
21	S481	PASS	523	10
22	S380	PASS	585	10
23	S151	PASS	667	10
24	S418	WITHDRAWN	-	WITHDRAWN
25	S339	PASS	527	10
26	S213	PASS	494	10
27	S167	PASS	747	10
28	S322	PASS	575	11
29	S449	PASS	631	11
30	S480	PASS	547	11
31	S323	PASS	663	11
32	S424	PASS	551	11
33	S476	WITHDRAWN	487	WITHDRAWN
34	T178	PASS	465	19
35	T171	PASS	409	19
36	T193	PASS	496	18
37	T177	PASS	417	18
38	T200	PASS	458	18
39	T342	PASS	462	18
40	T164	PASS	415	12
41	T208	PASS	456	18
42	T274	PASS	403	18
43	T185	PASS	469	18
44	T181	PASS	486	18
45	T273	PASS	457	4A
46	T445	PASS	433	17
47	T532	PASS	409	17
48	T188	PASS	449	17
49	T209	PASS	433	17

LOT	TAG	SEMEN TEST	WEIGHT (KG)	BULL PEN
50	T194	PASS	436	17
51	T525	PASS	427	17
52	T174	PASS	400	17
53	T569	PASS	432	17
54	T366	PASS	407	17
55	T575	PASS	432	15
56	T341	PASS	431	15
57	T330	PASS	460	15
58	T191	PASS	448	15
59	T187	PASS	419	15
60	T336	PASS	434	15
61	T192	PASS	421	15
62	T286	PASS	418	15
63	T207	PASS	389	16
64	T172	PASS	441	16
65	T240	PASS	427	16
66	T277	PASS	396	16
67	T201	QP	445	16
68	T244	PASS	377	16
69	T382	PASS	394	16
70	T180	PASS	443	13
71	T202	PASS	387	4A
72	T421	PASS	405	4A
73	T347	PASS	402	13
74	T280	PASS	430	13
75	T435	PASS	382	13
76	T196	PASS	431	13
77	T189	Pass	387	13
78	T175	PASS	376	13
79	T337	PASS	402	13
80	T494	PASS	380	13
81	T448	PASS	394	12
82	T275	PASS	407	12
83	T362	PASS	417	12
84	T210	PASS	393	12
85	T290	PASS	396	12
86	T523	WITHDRAWN	381	WITHDRAWN
87	T493	PASS	372	12
88	T426	PASS	383	12
89	T507	PASS	375	3
90	T247	PASS	378	3
91	T390	PASS	383	3
92	T170	PASS	372	3
93	T248	PASS	375	3
94	T167	PASS	375	3
95	T166	PASS	384	3
96	T451	PASS	382	3
97	T369	PASS	380	4
98	T361	PASS	373	4

LOT	TAG	SEMEN TEST	WEIGHT (KG)	BULL PEN
99	T449	PASS	396	4
100	T423	PASS	372	4
101	T199	PASS	407	4
102	T372	PASS	378	4
103	T335	PASS	414	4
104	T503	PASS	425	4
105	T437	PASS	378	5
106	T574	PASS	397	5
107	T357	PASS	403	5
108	T343	PASS	411	5
109	T349	WITHDRAWN	376	WITHDRAWN
110	T433	PASS	387	5
111	T436	PASS	370	5
112	T410	PASS	362	5
113	T227	PASS	355	5
114	T383	PASS	365	5
115	T235	PASS	373	8
116	T291	PASS	382	8
117	T246	PASS	371	8
118	T497	PASS	364	8
119	T242	PASS	377	8
120	T205	PASS	357	8
121	T243	PASS	368	8
122	T570	PASS	336	8
123	T380	PASS	393	7
124	T173	PASS	360	7
125	T443	PASS	367	1
126	T428	PASS	357	7
127	T491	PASS	372	7

LOT	TAG	SEMEN TEST	WEIGHT (KG)	BULL PEN
128	T438	PASS	370	7
129	T176	PASS	376	7
130	T442	PASS	369	1
131	T529	PASS	383	7
132	T425	PASS	372	7
133	T348	PASS	378	7
134	T233	PASS	370	6
135	T354	PASS	356	6
136	T272	PASS	387	6
137	T360	PASS	331	3
138	T510	PASS	355	6
139	T499	PASS	331	6
140	T281	PASS	352	6
141	T283	PASS	352	6
142	T508	PASS	365	1
143	T367	PASS	467	1
144	T195	PASS	440	1
145	T579	PASS	427	1
146	T184	PASS	421	1
147	T338	PASS	387	1
148	T197	QP	456	1
149	T271	QP	477	1
150	T522	PASS	357	1
151	T183	PASS	425	2
152	T490	QP	417	2
153	T206	WITHDRAWN	364	WITHDRAWN
154	T179	PASS	386	2
155	T339	PASS	351	2
156	T239	PASS	355	2

Semen Testing

All bulls have been semen tested by the industry leading team from Ced Wise Artificial Breeding Services.

Classifications of Semen Testing:

Pass (P)

Qualified Pass (QP)

A raw bull semen sample is classified as “pass” when motility is above 60% and normal forms greater than 70%. When samples were between 30 and 60% motility and 50–70% normal forms, they are classified as a “compensatory” (or qualified) pass (QP). The compensatory pass relies on there being the ability to have at least 10 million motile normal forms of spermatozoa in each straw to allow for conception. An outright failure is given to any sample with less than 30% progressive motility or 50% normal forms. This allows each ejaculate to be placed into a binary “pass” or “fail”.